

Genius Holding introduction

Innovative Applications in Sound and
Ultrasound Technology

Duncan Chen

dc@genius-gh.com

WhatsApp: +886926119246



Company Introduction

1

Genius

Genius Holding founded in 2016, provides creative ultrasound innovations to upgrade products with digitalized and wireless, and AI audio analysis for customers.



2

ezOxygen Smart Breath Checker-2019 CES Innovation Award

- ezOxygen personal spirometer with CES 2019 Innovation Award, Sport Fitness, and Biotech category.
- Strategic cooperation with **AstraZeneca** to launch “Health Lung App 呼吸肺健康”



3

Our Partnership

Provide Audio AI for smart manufacturing with **Chang Chun Petrochemical** (iEar Project), **Formosa Petrochemical corporation** (machine monitoring), **Asia Vital Components** (Product QA), etc.



Technology Overview

Passive Ultrasound Tech

Ultrasonic Data Transmission

AI Sound Analysis

1 Respiratory

Passive ultrasonic airflow determination: respiratory area (**ezOxygen**), anemometer, and gas flowmeter.

2 Smart Manufacturing

Audio Analysis with hearable and multiple sensors ultrasound wireless.

3 Fitness, VR/AR/MR, home security

Ultrasound Motion Tracking

4 Digitalized Medical Device

Ultrasound wireless medical devices, such as blood pressure, blood sugar, weight scale, etc

5 Non-Digitalized Medical Device

Ultrasound wireless medical devices, such as incentive spirometer, Intravenous (IV) Sets, etc.

Passive Ultrasound Technology




Airflow
transducer

Microphone

Algorithm &
display

Passive ultrasonic detection

The airflow generates a resonant frequency response through the metal whistle, and the signal to flow rate and volume are both calculated by the exclusive algorithm



Galton's whistle

Flow Direction

Ultrasonic Signal

Ultrasound airflow rate patent



Advantage



- Better accuracy than traditional method and no inertia effect or environment interference.
- **10,000 airflow rates** per second, v.s. traditional between **50~400 airflow rates** per second.)




Calibration-free.



Mouthpiece is easy to clean, and detachable




Ultrasonic Data Transmission




Ultrasonic Data Transmission

Sound wave communication- speaker and microphone.




Ultrasonic Data Transmission

Sound wave communication- speaker and microphone.



Ultrasonic Data Transmission


Sound wave communication- Position Trigger Module










Ultrasonic Data Transmission

Sound wave communication-Motion Detection


Real-time AI




Stop	Straight Punch		Uppercut		Lower Hook	
						

AI Sound Analysis


Applied to Equipment Maintenance and Finished Product QC Inspection



Sound Analysis -AI Validation Model



Time (2 sec)





Time (2 sec)

AI Model Specification :

- Classification Type : OK 、 NG
- OK Training Count : 410 pcs
- NG Training Count : 385 pcs
- Accuracy : 1.0
- Loss : 0.001719

Using audio data from verify to train the AI model to achieve testing.



Social Media/Customer Overview-Smart Factory

Digitimes News



https://www.digitimes.com.tw/tech/dt/n/shwnws.asp?CnId=13&cat=20&cat1=50&id=0000623141_BQ12VVDm3iV57S2IRUNRR&ct=d

Digitimes Media



【D Forum 2021智慧工廠論壇】傑尼斯：超聲波與聲音AI在智慧製造的應用

<https://www.youtube.com/watch?v=CGhu8rRvCkc>

Our Customer



<https://www.ccp.com.tw/>



<https://www.fpcc.com.tw/tw>



<https://www.avc.co/zh-tw/>

Comparison of Ultrasound Wireless Devices

AliveCor, Omron, and Genius

Products/year	Intend use	Differences	Common advantage
KardiaMobile (2014)	ECG, single lead	LTC6990 for voltage changed by fingers (18.5~19.5 kHz); distance <20-30 cm	Smartphone connect, App in iOS and Android, low battery consumption,
Omron MC-6800B (2021)	Body temperature	Frequency ranged between 18.5~19.5 kHz; IC(EPSON 17V14101 G2390234 IC + VRK272 IC); distance < 4 cm.	
Genius Ultrasound Incentive Spirometer (2023)	Lung function rehab	Customized crystal oscillators with specified kHz (18, 19.5, 20); distance <100~160 cm (by smartphones' sensitivity), with calibration to improve signals-noise ratio. Specified coding algorithm for other medical devices.	

Appendix - 1

AliveCor KardiaMobile ECG


The first & only ultrasonic ECG mobile device



- The original KardiaMobile. FDA-cleared, **single-lead** personal EKG that records medical grade EKGs right on your phone.
- The ultrasonic frequencies are ranged between 18.6~19.1 kHz.
- The device has a 10 mV peak-to-peak input dynamic range, a 300 samples per second sampling rate, and a 16-bit resolution.
- The device is powered by a 3V CR2016 coin cell battery that lasts for about 12 months of typical use.
- The device dimensions are 8.2 cm* 3.2cm*0.35 cm and it weighs 18 grams.
- The distance between the device and the smartphone depends on the placement of the device. If the device is attached to the back of the phone or the phone case, then the distance is minimal. If the device is not attached, then it can be placed within a few centimeters of the phone's microphone for optimal signal transmission.

Core ICs of KardiaMobile

LTC6990 & AD8502



LTC6990 is a **voltage-controlled oscillator (VCO) with a wide frequency range**, Key features of the LTC6990:

- Wide frequency range: 488 Hz to 2 MHz ; 4 frequency selection pins (FSEL) for 16 possible configurations
- Frequency adjustment via a single external resistor or voltage input ; CMOS output with 50% duty cycle
- Low supply voltage range: 2.25 V to 5.5 V ; Low supply current: 72 μ A at 100 kHz
- Temperature stability: ± 50 ppm/ $^{\circ}$ C (typical)

AD8502 is a precision, low power, **dual operational amplifier** (op-amp) require high accuracy and low power consumption, Key features of the AD8502:

- Low offset voltage: 1 μ V (typical) ; Offset voltage drift: 0.005 μ V/ $^{\circ}$ C (typical)
- Wide supply voltage range: 2.7 V to 5.5 V ; Low supply current: 1.1 mA per amplifier (typical)
- Rail-to-rail input and output ; Gain bandwidth product: 2 MHz (typical)
- Unity-gain stable

Omron音波通信体温計MC-6800B

2021 on Japan market ; 短時間で検温できる、約15秒の予測検温



- このことから、通信成功率が35%以上であれば、5秒以内に95%以上の確率で通信が完了すると考えられる。
- 表2は、各測定ポイント（奥行×3、高さ×2）における、通信成功率の分布をまとめたものである。目標とする通信成功率を達成できなかった機種は2~5%程度あるが、94%以上の機種で通信成功率80%以上という高い通信精度を達成している。
- この検証では、一度スマートフォンで採取したPCMデータをwavファイルに保存しておき、復調アルゴリズムをチューニング→wavファイルを用いて復調アルゴリズムの検証を繰り返すことで、復調アルゴリズムおよび各種パラメータのブラッシュアップにも寄与している。
- EPSON 17V14101 G2390234 IC + VRK272 IC

